What's Up in CT ?

BUC 07 Novembre 2024

Maxima Made in France

Scanner GE - Made in France

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Clinical Image gallery

GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.


*ECG-less Cardiac is 510(k) Pending at FDA. Not Available for Sale in the United States

Acute Chest Pain


Scan type	Axial 160 mm 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	60
Reconstruction	DLIR-H
kV	100
mA	503
CTDI	22.46
DLP, mGy-cm	359.33
mSv (*0.014)	5.03


Axial 0.625 mm



RCA

LCX

LAD

Multiphase Cine



Multiphase 4 Chamber



Multiphase Mitral Valve

Multiphase - Cardiomyopathy, Rule Out Coronary Artery Disease



Scan type	Axial 160 mm 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.28 sec
Estimated BPM	75
Reconstruction	DLIR-H
kV	100
mA	527
CTDI	29.59
DLP, mGy-cm	473.44
mSv (*0.014)	6.63


Multiphase Mitral Valve


Multiphase Pulmonary Valve

LAD


LCX


RCA

ECG-less cardiac acquisition estimated heart rate of 110 BPM

Scan type	Axial 160 mm 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	110
Reconstruction	DLIR-H
kV	100
mA	406
CTDI	9.87
DLP, mGy-cm	157.91
mSv (*0.014)	2.21


RCA

LAD





LCX

MIP LAD



TAVR

Scan type	Axial 160 mm 1 beat ECG-less cardiac (Full Cycle) Followed with Helical Aorta
Rotation time, s	0.28 sec
Estimated BPM	75
Reconstruction	DLIR-H
kV	100/100
mA	556/mod.
CTDI	20.99
DLP, mGy-cm	335.8
mSv (*0.014)	4.70


Post Stenting Follow Up

Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	60
Reconstruction	DLIR-H
kV	120
mA	315
CTDI	12.93
DLP, mGy-cm	206.89
mSv (*0.014)	2.90


LAD

LCX


RCA

MIP LAD

Preoperative AVR Evaluation


Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	51
Reconstruction	DLIR-H
kV	100
mA	625
CTDI	32.23
DLP	515.6
mSv (*0.014)	7.22


MIP


LCX


LAD


RCA

Aortic Valve w calcifications



Volume Illumination w TAVI measurements


Aortic Valve Sinus measurements

Full Beat Dynamic Capabilities


Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.28 sec
Estimated BPM	85
Reconstruction	DLIR-H
kV	120
mA	685
CTDI	34.11
DLP, mGy-cm	545.81
mSv (*0.014)	7.64


Multiphase Cine

Multiphase Aortic Valve


Multiphase Mitral Valve

RCA

LAD


LCX

MIP LAD

Calcified Aortic Valve


Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.28 sec
Estimated BPM	80
Reconstruction	DLIR-H
kV	100
mA	323
CTDI	11.54
DLP, mGy-cm	184.57
mSv (*0.014)	2.58

LCX

LAD



VR Aortic Valve


Multiphase Cine AV

CCTA – Myocardial Bridging of LAD

Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.28 sec
Estimated BPM	72
Reconstruction	DLIR-H
kV	100
mA	321
CTDI	12.52
DLP, mGy-cm	200.37
mSv (*0.014)	2.81


MIP LAD


LAD

RCA

Distal LAD Cross section


LCX

ECG-less cardiac acquisition estimated heart rate of 75 BPM


Scan type	Axial APEX
	1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	75
Reconstruction	DLIR-H
kV	100
mA	624
CTDI	22.92
DLP, mGy-cm	366
mSv (*0.014)	5.1

MIP LAD


LAD


RCA

Evaluation Of Heart Failure Etiology


Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	60
Reconstruction	DLIR-H
kV	100
mA	937
CTDI	38.91
DLP, mGy-cm	622.58
mSv (*0.014)	8.72


Aortic Valve



Aortic Valve

RCA

LAD



MIP LAD


Mitral Valve

55 BPM


Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	55
Reconstruction	DLIR-H
kV	100
mA	458
CTDI	20.46
DLP, mGy-cm	286.49
mSv (*0.014)	4.01


RCA


LMB


RIB

LAD


LCX

Axial 0.625mm

ECG-less cardiac acquisition estimated heart rate of 60 BPM


Scan type	Axial APEX
	1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	60
Reconstruction	DLIR-H
kV	100
mA	527
CTDI	23
DLP, mGy-cm	376
mSv (*0.014)	5.3

RCA

LCX

LAD

Angina Pectoris At High Risk

Scan type

Rotation time, s

Estimated BPM

Reconstruction

kV

mA

CTDI

DLP, mGy-cm

mSv (*0.014)

Axial APEX

1 beat ECG-less cardiac (Full Cycle)

0.28 sec

75

DLIR-H


100

863

44.25



708

9.91


Tortuous Coronary Arteries scanned with estimated BPM of 75


Scan type	Axial APEX
	1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.28 sec
Estimated BPM	75
Reconstruction	DLIR-H
kV	120
mA	624
CTDI	48.77
DLP, mGy-cm	780.32
mSv (*0.014)	10.92


LAD

LCX

LMB


RCA


Lumen RCA

History Of Multivessel Coronary Artery Disease. Surveillance

Scan type	Axial APEX 1 beat ECG-less cardiac (Full Cycle)
Rotation time, s	0.23 sec
Estimated BPM	65
Reconstruction	DLIR-H
kV	120
mA	775
CTDI	51.97
DLP, mGy-cm	831.49
mSv (*0.014)	11.64

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



Ascend



Improved contrast enhancement for PE with True Enhance DL*

Scan type Helical 40 mm
Rotation time, s 0.5
Pitch 1.375
Slice, mm 1,25
Reconstruction AR50
kV 120
WW WL 800 150

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




AWS IA et Deep Learning

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
GE HealthCare

GE HealthCare Photon Counting CT Vision

Build a general-purpose, clinical PCCT to combine breakthroughs in Spectral and Spatial with Coverage, Speed, Low Dose and Workflow.

At the same time!


Deep silicon in perspective

B	C	N
Boron	Carbon	Nitrogen
10.81	12.01	14.01
Al	Si	P
Aluminum	Silicon	Phosphorus
26.98	28.09	30.97
Ga	Ge	As

Why Silicon?


- Purest material and perfect crystal
- 40 years ahead of any other semiconductor
- Spectral sensitivity without compromising spatial resolution
- Futureproof
- Continuous and rapid development of available processes
- Environmentally safe
- Cd banned in the EU: Restriction of Hazardous Substances Directive 2006, exemption only given 7 years at a time, when no alternative exist.

The purity of Si is 100,000 times higher than CdTe.^{1,2}

Semiconductor capacity expansion needs trillions of dollars in investment - GlobalFoundries CEO

Jun. 23, 2022 11:07 AM ET GLOBALFOUNDRIES Inc. (GFS) AMD, INTC, NVDA | By: Brian Stewart, SA News Editor 65 Comments

Phonlamai/PhotoStock via Getty Images

GlobalFoundries CEO Thomas Caulfield said Thursday that the semiconductor sector will require "trillions of dollars" of investment to reach a point where it can meet global demand.

PCCT – Deep Silicon Detector


“Having developed the world’s first photon counting CT prototype in 2006, GE came to understand the strengths and limitations of the cadmium-based detectors that have been the primary focus of photon counting CT research, and upon which that prototype relied, says Jonathan S. Maltz, GE’s Chief Scientist for CT.

“As well as high spatial resolution, the enhanced spectral resolution realized by Deep Silicon technology has the potential to allow detection of subtle changes in tissue properties, expanding the potential of CT images,” says Maltz.

[Read Less](#) ^


Table 1. Comparison of properties of cadmium-based and silicon-based detectors. Compared to a CdTe/CZT detector, constructing a full-area silicon detector system is more complex because of the edge-on sensor geometry. FWHM, full width at half maximum.

Detector material	CdTe/CZT	Silicon
Preferred orientation	Face-on	Edge-on
System construction complexity	Less complex	More complex
Typical lowest threshold	20–25 keV	5–10 keV
Typical energy resolution (FWHM)	5–10 keV	3.5–5.4 keV
K-fluorescence escape probability	High	Low
Compton-scatter probability	Low	High
Mobility-lifetime product	Low	High
Amount of polarization at high fluence rate	High	Low
Material processing technology	Under development	Mature


Deep Silicon detector innovations to realize Spectral+Spatial breakthroughs

Edge-on geometry with multiple 'counters' to count each x-ray photon separately


Deep Silicon technologies


Multi-Energy-Bin for better energy separation and high spectral resolution

Conventional technologies

Probability distribution of original photon energy: CdTe

PCCT full clinical potential requires profound improvement in Spectral + Spatial and the combination of both¹

- Earlier Cancer Diagnosis
- Opportunistic Screening
- Increased Anatomic Delineation
- Superior Quantification
- Lesion Characterization

Oncology


- Superior grey white differentiation
- Increased vascular assessment with artifact reduction
- Increased CNR for small lesion detection
- Enhanced Lesion Characterization

Neuro


Cardiac

- Superior plaque delineation
- Increased accuracy of stenosis quantification
- Minimization of Calcium Blooming
- Morphological Plaque characterization
- Myocardial assessment

Spectral
+
Spatial

MSK

- Increase Anatomic Delineation
- Disease Characterization
- Accurate Crystal Separation
- Quantitative Imaging

1. Experts' consensus of GE PCCT Medical Advisory Board Steering Committee.


Generations Deep Silicon

Gen 1 Prototype 2016

Gen 2 Prototype 2021

Gen 3 Prototype 2023

Gen 4 ME System

Proof of concept

Technology Demonstration

Enhanced Technology Demonstration

Final product, serviceable system

for investigational
in-vitro/vivo
studies

for investigational
studies on
human subjects

for investigational
studies on
human subjects

Full capacity
clinical system

Ultra-high-resolution

Volumetric Hi-Res spectral

High pitch helical

ECG-gated cardiac

Routine clinical

PCCT – Deep Silicon Detector


Figure 4. A comparison of the high-resolution imaging performance on the 28 lp mm⁻¹ (A), 22 lp mm⁻¹ (B) and 14 lp mm⁻¹ (C) patterns of a CATPHAN 700 phantom (CTP714, The Phantom Laboratory) for (a) a conventional (single-energy) CT system with 0.7 mm focal spot, (b) a prototype silicon-based PCD with 0.4 mm focal spot. Both images were acquired with axial mode at 120 kVp with comparable dose and reconstructed with a high-resolution kernel. Adapted with permission from da Silva *et al* (2019). © 2019 Society of Photo-Optical Instrumentation Engineers (SPIE).


Figure 1. An image of an excised human heart with a calcified coronary artery containing iodinated gd. The leftmost image is a 67 keV virtual monoenergetic image. Two-material basis decomposition has been performed with the basis pairs: calcium (water), iodine (water), calcium (iodine), and iodine (calcium). The image set was acquired using a state-of-the-art dual-energy CT scanner. The slice thickness was 0.625 mm; the focal spot was 1.2 mm; and the x-ray tube was operated at 120 kVp and 384 mAs. The images were acquired at the Karolinska Hospital, Stockholm, Sweden. Adapted from Grünberg *et al* (2020) (© 2020 Springer Nature Switzerland AG. Part of Springer Nature.). Adapted with permission of Springer. CC BY 4.0.


Figure 2. An image of an excised human heart imaged with a prototype silicon-based spectral PCD. The leftmost image is a 67 keV virtual monoenergetic image. Three-material basis decomposition was performed with water, calcium, and iodine as bases. An overlay image was formed in which the regions of calcium and iodine are colored red and green respectively. The image was acquired at 120 kVp using a focal spot 0.4 mm, and reconstructed with a 0.625 mm slice thickness. The dose was matched to that of the dual-energy image in figure 1. Adapted from Grünberg *et al* (2020) (© 2020 Springer Nature Switzerland AG. Part of Springer Nature.). Adapted with permission of Springer. CC BY 4.0.

Initial Clinical Images From a Second-Generation Prototype Silicon-Based Photon-Counting Computed Tomography System

Hakan Almqvist MD,^{a,b} Dominic Crotty PhD,^c , Sven Nyren MD, PhD,^{d,e} Jimmy Yu MD,^{d,e} Fabian Arnberg-Sandor MD, PhD,^{a,f} Torkel Brismar MD, PhD,^{g,h} Cedric Tovattⁱ, Hugo Linder MS^j, Jose Dagotto BS^l, Erik Fredenberg PhD^l, Moa Yveborg Tamm PhD^l, Paul Deak PhD^l, Michail Fanariotis MD, PhD,^k Robert Bujila PhD^l, Staffan Holmin MD, PhD,^{a,f}

[Show more](#) 